



MODEL M1204

Capacity:	Holds four bronchoscopes, Cystoscopes, or Uretheroscopes	Separator Rack:	Center of cabinet Two slots for each scope One slot for insertion tube-prevent cross contamination One slot for support of connector One-piece 3/4" polyethylene construction Cam-lock fasteners secure to cabinet side and back walls
Ext. Dim:	13"D x 21.6"W x 77"H	Back Wall:	Padded vinyl to prevent damage to distal tips of scopes Inset 1/2" to permit ambient air circulation
Weight:	130 lbs.	Drip Tray:	Removable polystyrene drip tray
Color:	Architectural White	Air Vents:	Top and bottom of cabinet allow ambient air circulation
Construction:	3/4" Melamine Back Wall 5/8" Melamine		
Trim:	Brushed aluminum T-Mold and Kickplate		
Door:	FlexLok™ Tambour door rolls up around counter balance key locked		
Cabinet Top:	Powder coated steel - totally slanted to prevent dust collection		
Bottom Surface:	Finished to prevent moisture absorption		
Scope Rack:	One piece 3/4" polyethylene construction Polyethylene brackets provide added support for weight of scopes Cam-lock fasteners secure rack to cabinet side and back walls Closed cell foam cushions protect scope heads Holds four small diameter scopes		



MASS Medical Storage

10035 Lakeview Ave. Lenexa, KS 66219

913-438-8835 Fax 913-438-2275

www.massmedicalstorage.com

M1204