



MODEL M2000HG

Capacity: Eight flexible endoscopes

Ext. Dim: 21.6"D x 41"W x 93.15"H

Int. Dim: 20"D x 39.37"W x 88.65"H

Weight: 275 lbs.

Color: Architectural White

Construction: 3/4" Melamine
Back wall 5/8" melamine

Trim: Brushed aluminum T-Mold and Kickplate

Door: Hinged doors with 1/8" tempered glass inserts on each side
Optional: Solid hinged doors on both sides
Hinged door w/ 1/8" tempered glass insert on one side - solid hinged door on other side
Locking doors on one or both sides

Cabinet Top: Flat top for installation in wall cutout

Bottom Surface: Finished to prevent moisture absorption

Scope Rack: One piece 3/4" polyethylene construction
Polyethylene brackets provide added support for weight of scopes
Cam-lock fasteners secure rack to cabinet side and back walls
Closed cell foam cushions protect scope heads
Holds eight flexible endoscopes

Insertion Tube Separator Rack: Center of cabinet - prevent cross-contamination
One piece 3/4" polyethylene construction
Cam-lock fasteners secure to cabinet side and back walls

Back Wall: Padded vinyl to prevent damage to distal tips of scopes

Drip Tray: Removable polystyrene drip tray

Air Vents: Bottom of cabinet - front and back; filtered
Allow ambient air circulation



MASS Medical Storage

10035 Lakeview Ave. Lenexa, KS 66219

913-438-8835 Fax 913-438-2275

www.massmedicalstorage.com

M2000HG