

**Model Numbers:** SAS-9M/P/S-LHG, SAS-18M/P/S-LHG

**Color:** Bright White Melamine  
Bright White Powder Coated Metal  
Stainless Steel

**Cabinet Filter Efficiency:** Removes 99.97% of Particles (Up to 0.3µm)

**Air Flow Rate Filter Efficiency:** 481 L/min

**Capacity:** 9 or 18 Flexible Endoscopes

**Ext. Dimensions:** 18 3/4"D X 40 7/8"W X 88 5/8"H

**Construction:** Melamine, Powder Coated Metal, or Stainless Steel  
Locking Glass Doors  
Padded Lower Walls (Distal Tip Protection)

**Air Supply:** HEPA Filter Fan (Cabinet)

**Electrical Supply:** 120 V Wall Plug-3 A (10' Cord Exiting Top Right of Cabinet)

**Options:** Door Activated LED Lighting  
Keyless Lock (Mechanical or Digital)  
Temperature and Humidity Gauge

## Features



Measurements

120V Hospital Grade Plug  
Exits Top Right Side. Approx. 10' Length



