

MASS

Reduction
of Nosocomial Infection

Protection
of Valuable Medical Equipment

Progression
of Clinician Safety

MASS[™] truAIR[™]

HEPA-Grade Forced Air Endoscope Drying Cabinets

The **NEW** MASS[™] truAIR[™] Endoscope Drying Cabinets provide industry leading storage and drying.

Offering HEPA-filtered, airflow directly into the cabinet itself and continuous positive air pressure through each endoscope lumen at 99.97% filtration up to 0.3µm

Our expanded line of drying cabinets make compliance with AORN, SGNA, AAMI, and The Joint Commission accessible to more facilities than ever before.

FEATURES | MASS[™] truAIR[™]


- High Flow Rate Channel Drying with 6LPM per Scope
- 45dB Volume Equivalent to Light Rain
- Self-Contained Air System
- 2 Years Maintenance-Free Continuous Use Air Compressor
- Rotating Scope Hooks
- Reusable Drying Connectors including Auxiliary Water
- Antimicrobial Powder Coated Metal
- 4 or 8 Scope Capacity Half Width
93 3/4" H x 21 5/8" W 18 3/4" D
- 9 or 18 Scope Capacity Full Width
93 3/4" H x 40 7/8" W 18 3/4" D




ALSO AVAILABLE
GI Transport Carts
and Accessories




1 CLICK



2 CONNECT



3 DRY



Call 800.593.1900 for a Quote

MASSMedicalStorage.com