



**Reduction**  
of Nosocomial Infection

**Protection**  
of Valuable Medical Equipment

**Progression**  
of Clinician Safety

# MASS<sup>™</sup> truAIR<sup>™</sup>

Hospital Air Endoscope Drying Cabinets

The **NEW** MASS<sup>™</sup> truAIR<sup>™</sup> Endoscope Drying Cabinets provide industry leading storage and drying.

Offering HEPA-filtered, pressurized airflow directly into the cabinet itself and continuous positive air pressure through each endoscope lumen at 99.97% filtration up to 0.01µm

Our expanded line of drying cabinets make compliance with AORN, SGNA, AAMI, and The Joint Commission accessible to more facilities than ever before.

## FEATURES | MASS<sup>™</sup> truAIR<sup>™</sup>


- Connects to Hospital Air System
- Rotating Scope Hooks
- Reusable Drying Connectors including Auxiliary Water
- Air Manifold with Quick Disconnect
- Locking Glass Doors
- Center Separator Scope Rack
- Removable Drip Tray
- Antimicrobial Powder Coated Metal
- 4 or 8 Scope Capacity Half Width  
93 3/4" H x 21 5/8" W 18 3/4" D
- 9 or 18 Scope Capacity Full Width  
93 3/4" H x 40 7/8" W 18 3/4" D




**ALSO AVAILABLE**  
GI Transport Carts  
and Accessories




**1** CLICK



**2** CONNECT



**3** DRY



Call 800.593.1900 for a Quote

MASSMedicalStorage.com